

## TOP CHECK 2 - توب تشيك ٢

**TOP CHECK 2** is a semi-automated photo-optical coagulation and turbidity measurement system. It's designed to perform clotting and turbidimetric assays.

It has auto-start function in all programs initiated with sample addition. All advantages are put in a unique maintenance free design.



### **Principle of operation :**

When specimen is put into optic well, the optical system measures its optical absorbance and the rate of change during the test interval, then data is processed by the computation unit with a presaved routine to generate the results that are finally shown on LCD.

### **Measurement Methods :**

- 1- Coagulation ( Time % INR Ratio Conc.)
- 2- Single standard fixed rate (  $\Delta$ Abs. Conc.)
- 3- Multi-standard fixed rate (  $\Delta$ Abs. Conc.)

### **Measureable Tests :**

- PT PTT TT FIB. Factors Lupus A.
- D.Dimer Ferritin HbA1c CRP ASO RF.
- Syst.C IgA IgE IgG IgM C3 C4 Vit.D.
- C.Trop.  $\mu$ ALB. Lipo P. C.Plasmin Haptoglobin G.

### **Software :**

- Easy test selection (TEST Btn)
- Easy calibration (PG, Up, Dn, EnterBtns)
- Auto Start function (RESET Btn)
- Timer and manual Start (START Btn)
- Useful test summary help is shown on LCD second line, eg. [ 50u Rgt + 25u S ]

### **Incubation :**

All wells are electronically thermostated to 37.1°C.

- 6 wells for samples incubation.
- 2 wells for reagent / tube incubation.

### **Minimum sample volume :**

- 75 $\mu$ L (25+50) minimum test volume.

### Test List (TOP CHECK 2 - ver 2.12)

No.	Test Name	Reagents & Procedure						T.min	T.max	Test Method	Default values			
		Rgnt	Smpl	R1	R2	Sample	R1	Sample	R2	Delay	Fixed.T	Normal	STD.	Factor
1	PT	50	25							6	300	12.9	ISI=1.5	
2	PTT						25	25	25	6	300	32		
3	TT	50	50							6	300	15		
4	FIB.	25	50							6	300	200-400	MTSD	
5	FAC	50	25							6	300		MTSD	
6	Lupus A.						25	25	25	6	300		42	
7	D.Dimer						45	10	15	40	150	< 0.5	3.65 /2	10.5
8	Ferritin			80	20	10				20	150		525	1600
9	HbA1c						75	5	25	20	150	< 80	265	750
10	CRP			100	25	5				20	120	< 6	63 /2	100
11	ASO			100	25	5				20	120			1000
12	Rhum.F.			100	25	5				20	150			250
13	Syst.C			100	20	5				20	150			16
14	IgA	100	5							20	150			1000
15	IgE			75	25	10				20	150			500
16	IgG	100	5							20	150			4000
17	IgM	100	10							20	150			750
18	C3	100	5							20	150			600
19	C4	100	5							20	150			160
20	Vit.D						80	5	20	20	150	20-40	50	160
21	C.Troponin						75	10	25	20	150			50
22	Micro ALB.						80	5	15	20	150			400
23	Lipo P.			100	15	5				20	150			1600
24	C.Plsmn						100	5	15	20	150			160
25	Hapto G.						100	5	15	20	150			800
26	FIB. FxT						100	5	15	20	150	200-400	300	100
27	Multi.STD 1									20	240	FxT.	MTSD	
28	Multi.STD 2									20	240	3 points	MTSD	
29	Multi.STD 3									20	300	Multi	MTSD	
30	Multi.STD 4									20	300	STD.	MTSD	