












Chem-Tech



www.chemtech-eg.com





## oven

CHEM.TECH Dry heat oven is used in many applications. Generally, the temperature range is from 50 to 300 degrees Celsius. There is a temperature control thermostat.

Process applications for laboratory ovens can be for drying, Polyimide baking, sterilizing, and other industrial laboratory functions.

## **Features**

- 1- LCD digital display enable users to monitor the chamber temperature at any given time.
- 2- Ability to set operating time of heater.
- 3- Temperature range from ambient temperature to 300°C.
- 4- Safety: Visual and auditory alarm indicator alerts users of temperature abnormal conditions.
- 5- External: coated with electrostatic.
- 6- Internal: high quality stainless steel.
- 7- removable Racks.
- 8-Easy cleaning and optimal airflow.
- 9-Keyed on/off switch.



Chem Tech Oven			
Model	OCT.25	OCT.40	OCT.60
capacity	25 L	40 L	60 L
Max Temp	300° c	300° c	300° c
Heater	400 w	600 w	800 w
Dimentions	53*44*51	53*51*56	53*60*62
voltage		220 v	