

FILIPS 400

Semi- Auto Biochemistry Analyzer



Technical Parameters

Light Source : Halogen Lamp 6 V, 10 W

Screen Touch , 7inch , Keyboard , Mouse ,USB

Suction System & Cuvette System. Injection: 0~6000ul

Interface : Standard Rs-232 interface , 4 USB interface,

FILIPS 400

Semi-Auto Biochemistry Analyzer



Introduction

Multiple measurement and multiple calculation methods

Store 500 test items and 10000 test results

Each item can be set for two quality-control

Self-check function for alarming to light route etc

With standard interface for communication

convenient for data transmitting and processing

The Semi-auto Biochemistry Analyzer Measure

Biochemical Index by analyzing

body fluid, then combines With other clinical

information, to help diagnose disease, evaluate,

organs function, identify disease gene

determine the norm for future therapy.

